

# INDIPRO

## Rapid Test for determination of protein residues

---

### Contents:

- 1 aluminium container with 60 test strips
- 1 bottle INDIPRO-1
- 1 bottle INDIPRO-2

INDIPRO is a rapid test for determination of protein residues on working surfaces and tools having contact with foodstuffs.

### Safety precautions:

Information regarding safety can be found on the box' label and in the safety data sheet. You can download the SDS from [www.mn-net.com/SDS](http://www.mn-net.com/SDS).

### General indications:

Remove only as many test strips as are required. Close the container immediately after removing a strip. Do not touch the test field.

### Instructions for use:

1. Add 1 drop of reagent INDIPRO-1 on the test field of the test strip.
  2. Draw test strip with the moistened test field approximately 10 cm under low pressure over the surface to be tested.
  3. Add 1 drop of reagent INDIPRO-2 on the test field.
- A colour change from yellow to green indicates the presence of protein residues. The surface has to be cleaned.

### Limit of sensitivity:

50 µg BSA (bovine serum albumin) gives a clear green colouration.

### Interference:

Residues of desinfectants can cause also a green colouration of the test field.

### Disposal:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from [www.mn-net.com/SDS](http://www.mn-net.com/SDS).

### Storage:

Avoid exposing the strips to sunlight and moisture. Keep container cool and dry (storage temperature 4–30 °C).

If correctly stored, the test kit may be used until the use-by-date printed on the packaging.

**CTL SCIENTIFIC SUPPLY CORP.** 1016-3 Grand Boulevard, Deer Park, NY 11729

Tel: **631-242-4249**

Web: [www.ctlscientific.com](http://www.ctlscientific.com)

Manufacturer: **Macherey-Nagel GmbH & CO. KG**

Rev: 2025-03