

Manganese

**Test kit for performing colorimetric tests
on manganese ions in surface water and sewage**

Method:

Determination of total manganese with formaldoxime

Measurement range:

0.1–1.5 mg/L Mn

Contents of test kit (*refill pack):

sufficient for 70 tests

30 mL Mn-1*

30 mL Mn-2*

4 g Mn-3*

1 measuring spoon 70 mm*

2 screw-plug measuring glasses

1 slide comparator

1 color chart

1 plastic syringe 5 mL

1 instructions for use*

Hazard warning:

Information regarding safety can be found on the box' label and in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

Instructions for use:**a) colorimetric determination with color chart**

also refer to the pictogram on the back of the color chart

1. Pour a **5 mL water sample** (the pH value of the sample must be between pH 1 and 13) into each of the measuring glasses using the plastic syringe.

Place a measuring glass on position A in the comparator.

Only add the reagent to measuring glass B.

2. Add **5 drops of Mn-1**, seal the glass and mix.

3. Add **7 drops of Mn-2**, seal the glass and mix. Wait **1 min**.

4. Add **1 level measuring spoonful of Mn-3**, seal the glass and shake well.

5. Open the glass after **5 min** and place it on position B in the comparator.

6. Slide the comparator until the colors match in the inspection hole on top. Check the measurement reading in the recess on the comparator reed. Mid-values can be estimated.

7. After use, rinse out both measuring glasses thoroughly and seal them.

b) photometric determination

The reagents are also suitable for **photometric evaluation**. Please refer to the separate instructions for photometric performance.

The method can be applied also for the analysis of sea water.

Disposing of the samples:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

Interferences:

The following ions will not interfere:

≤ 0,1 mg/L Ni²⁺; ≤ 2 mg/L Cr, Co²⁺; ≤ 20 mg/L Cu²⁺, Fe³⁺;

≤ 50 mg/L Al³⁺, PO₄³⁻

Storage:

Store the test kit in a cool (< 25 °C) and dry place.