

REF 985047

en

Test 0-47

03.23

NANOCOLOR® Nonionic surfactants 15

Method:

Nonionic surfactants (ethoxylates with 3 to 20 ether bridges) react with an indicator (TBPE) to form a complex, which is extracted with dichloromethane.

Range:	0.3 – 15.0 mg/L Triton® X-100
Wavelength (HW = 5 – 12 nm):	620 or 610 nm
Reaction time:	2 min (120 s)
Reaction temperature:	20 – 25 °C

Contents of reagent set:

20 test tubes Nonionic surfactants 15

Hazard warning:

Test tubes contain dichloromethane 90–100 % and methanol 3–10 %.

H341, H351 Suspected of causing genetic defects. Suspected of causing cancer.

P201, P280 sh Obtain special instructions before use. Wear protective gloves/eye protection.

For further information ask for a safety data sheet.

Interferences:

Cationic surfactants cause high results.

Anionic surfactants cause low results, depending on the kind of the anionic surfactant.

The following ions will not interfere: < 1000 mg/L K⁺, Na⁺, Cl⁻; < 500 mg/L NH₄⁺, NO₃⁻, SO₄²⁻;

< 200 mg/L Ca²⁺, Mg²⁺, Ni²⁺, Zn²⁺; < 100 mg/L Fe³⁺, Al³⁺; < 50 mg/L Cu²⁺; < 10 mg/L NO₂⁻.

The method cannot be applied for the analysis of sea water.

Procedure:

Requisite accessories: piston pipette with tips

Open test tube, add

4.0 mL test sample (*the pH value of the sample must be between pH 4 and 9*), close and **shake gently for 2 min** (*with a frequency of 2–3 times per second*).

Clean outside of test tube and measure after 2 min (*wait for phase separation*).

Measurement:

For MACHEREY-NAGEL photometers see manual, test 0-47.

Note:

The calibration curve is calculated as Triton® X-100. For measuring other nonionic surfactants verify calibration curve by measuring standard solutions.

Photometers of other manufacturers:

For other photometers check whether measurement of round glass tubes is possible. Verify calibration curve for each type of instrument by measuring standard solutions.

MACHEREY-NAGEL GmbH & Co. KG · Valenciener Str. 11 · 52355 Düren · Germany

Tel.: +49 24 21 969-0 · info@mn-net.com · www.mn-net.com

PD 14122 / A031351 / 985047 / 033xx